



Solidification Point of Benzene



ASTM D852

Solidification Point of Benzene

This test method covers the determination of the solidification point of benzene.

Art. LT/SP-237100/M

Solidification Point According to ASTM D852

- Dewar made in glass
- Jacket tube
- Test tube (Benzene container)
- Cover for jacket tube with opening for dry ice
- Stopper for test tube
- Stopper for thermometer and stirrer
- Stirrer: stainless steel wire

Accessories

- T-AS112C: thermometer ASTM 112 C +4°C ... +6°C

Spare Parts

- LAB-182-341: dewar
- LAB-152-342: jacket tube 25 × 150 mm
- LAB-152-343: test tube 15 × 125 mm, pack 10 pcs.
- LAB-152-344: cover for jacket tube
- LAB-152-345: stopper for test tube
- LAB-152-346: stopper for thermometer
- LAB-152-347: stainless steel wire