



Wax Melting Point



ASTM D87
IP 55

Melting Point of Petroleum Wax (Cooling Curve).

This test method covers the determination of the melting point (cooling curve) of petroleum wax. It is unsuitable for waxes of the petrolatum group, microcrystalline waxes, or blends of such waxes with paraffin wax or scale wax.

Art. LT/WM-209000/M Wax Melting Point ASTM D87

- Air bath with brass well
- Glass test tube $\varnothing 25 \times 100$ mm calibrated to 50.8 mm
- Cork stoppers
- Water bath

Accessories

- LAB-100-332: digital stopwatch
- T-AS14C: thermometer ASTM 14C - IP 17C
- T-AS14F: thermometer ASTM 14F - IP 17F
- T-AS34C: thermometer ASTM 34C - IP 21C
- T-AS34F: thermometer ASTM 34F

Spare Parts

- LAB-102-091: calibrated dish, pack of 10
- LAB-102-092: cork